

MERCK

Having trouble with electrophoresis? Merck has the solution!

mPAGE® Bis-Tris precast gels

mPAGE® Bis-Tris precast gels provide excellent protein resolution at a budget-friendly price.



mPAGE® LUX CASTING SYSTEM

Cure a mini gel in 90 seconds with the mPAGE® Lux Casting System compared to conventional methods that take more than 90 minutes.

Our quick casting workflow with UV irradiation expedites curing of mPAGE® Lux Bis-Tris reagents while Bis-Tris gel chemistry allows for shorter run times compared to conventional methods.



**01ST Center Bldg. B503
13:00 – 17:00
23 October 2024**

Merck Ltd.

Phone: 03-4531-1145

Email: mj_is@merckgroup.com

Web:

<https://www.merckmillipore.com/JP/en>

<https://www.sigmaaldrich.com/JP/en>