

Demonstration at OIST

Confocal Laser Scanning Microscope

FLUOVIEW™ FV5000

**Just
Launched!**

SilVIR Next-Generation Detector

2K resonant & 8K galvo scanning

Automatic correction collar

Simple, intuitive Software

AI-assisted workflows



Product Information

Mouse brain slice expressing 7-color Tetbow cleared with SeeDB2 (in utero electroporation).
Sample courtesy of: Drs. Satoshi Fujimoto and Takeshi Imai, Graduate School of Medical
Sciences, Kyushu University.

Date & Time 4-6Feb 2026

4 Feb (Wed) , 5Feb (Thu) 10:00-17:30

6 Feb (Fri) 10:00-12:00

Location: OIST Lab1C352

“Scan to Register”



Contact : Tamaki(EVIDENT)masashi.tamaki@evidentscientific.com

EVIDENT